

# KERN 317-090-400



E1, E2 Plastic box for individual weights 500 g



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Plastic box

## Measuring System

Case suitable for OIML classes	E1, E2
Case suitable for nominal weight values	500 g
Suitable weight shape for case	Cylindrical/compact

## Construction

Material	plastic
Dimension (Ø×H)	61×94 mm

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	61×61×94 mm
Net weight	0,02 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

## Pictograms

STANDARD

