

KERN 317-120-100



E1+E2+F1 Wooden box for individual weights, lined 2 kg



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	E1+E2+F1
Case suitable for nominal weight values	2 kg
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	85×85×123 mm
Material	beech

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	85×85×123 mm
Net weight	0,35 kg
Shipping method	Parcel service
Net weight approx.	0,35 kg
Gross weight approx.	0,35 kg
Shipping weight	0,35 kg

Pictograms

STANDARD

