

# KERN 317-140-100



E1+E2+F1 Wooden box for individual weights, lined 10 kg



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

## Measuring System

Case suitable for OIML classes	E1+E2+F1
Case suitable for nominal weight values	10 kg
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (W×D×H)	126×126×216 mm
Material	beech

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	126×126×216 mm
Net weight	0,9 kg
Shipping method	Parcel service
Net weight approx.	0,90 kg
Gross weight approx.	0,90 kg
Shipping weight	0,9 kg

## Pictograms

STANDARD

