

KERN OZB-A5603



Planachromatic objective lens 1 x



Category

Brand	Optics
Product category	Optical Instrument
Product group	Microscope objective lens

Construction

Anti fungus	✓
-------------	---

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	55×55×25 mm
Net weight	0,5 kg
Shipping method	Parcel service
Net weight approx.	0,50 kg
Gross weight approx.	0,70 kg
Shipping weight	0,7 kg

Objective

Objective quality	Plan achromatic
Objective magnification	1 x

Pictograms

STANDARD

