

# KERN 303-054

Weight set, E1, 1 mg - 500 g, Knob, stainless steel polished, in plastic case



## Optional cases

Case - Plastic article number 313-052-400

## Product information

GTIN/EAN number 4045761291607  
REACH classification no prohibited substances

## Pictograms

### STANDARD



### OPTION



## Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	303

## Measuring System

Weight OIML class	E1
Test weight - nominal value	1 mg - 500 g
OIML-compliant	✓

## Construction

Dimension (W×D×H)	250×210×75 mm
Shape weight	Knob
Material	stainless steel polished
Number of weights	24
Case of weight set	in plastic case
Adjustment method	Polishing

## Services (optional)

Article number for DAkks calibration	963-204
--------------------------------------	---------

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	250×210×75 mm
Shipping method	Parcel service
Net weight approx.	3,2 kg
Gross weight approx.	3,2 kg
Shipping weight	3,19 kg