

KERN 304-036

KERN

Weight set, E1, 1 g - 100 g, Knob, Stainless steel polished, in aluminium case



Optional cases

Case - Aluminium article number 314-030-600

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	304

Measuring System

Test weight - nominal value	1 g - 100 g
OIML-compliant	✓
Weight OIML class	E1

Construction

Dimension (W×D×H)	165×125×75 mm
Shape weight	Knob
Material	Stainless steel polished
Case of weight set	in aluminium case

Services (optional)

Article number for DAkkS calibration	963-216
Article number for guaranteed plus-deviation	921-216

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	165×125×75 mm
Shipping method	Parcel service
Net weight approx.	0,65 kg
Gross weight approx.	0,90 kg
Shipping weight	0,9 kg