

KERN 304-06

KERN

Weight set, E1, 1 g - 1 kg, Knob, stainless steel polished, in wooden case



Optional cases

Case - Wood article number 313-060-100

Product information

GTIN/EAN number 4045761002586
REACH classification no prohibited substances

Pictograms

STANDARD



OPTION



Category

Brand	KERN
Product category	Weight set
Product group	OIML E1
Product family	304

Measuring System

Weight OIML class	E1
Test weight - nominal value	1 g - 1 kg
OIML-compliant	✓

Construction

Dimension (W×D×H)	230×180×110 mm
Shape weight	Knob
Material	stainless steel polished
Number of weights	13
Case of weight set	in wooden case
Adjustment method	Polishing

Services (optional)

Article number for DAkkS calibration	963-219
Article number for guaranteed plus-deviation	921-219

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	230×180×110 mm
Shipping method	Parcel service
Net weight approx.	4,4 kg
Gross weight approx.	4,4 kg
Shipping weight	4,21 kg