

# KERN 308-39

# KERN

Milligram weights, E1, 500 mg, Wire, stainless steel polished



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML E1
Product family	308

## Measuring System

Conformity assessment possible	✓
Test weight - volume (cm <sup>3</sup> )	0,0625
Test weight - volume uncertainty (cm <sup>3</sup> )	0,0013
Test weight - Density (kg/m <sup>3</sup> )	8000
Test weight - Density uncertainty (kg/m <sup>3</sup> )	160
Test weight - Density [Min] (OIML)	6300
Test weight - Density [Max] (OIML)	10900
Test weight - nominal value	500 mg
OIML-compliant	✓
Weight OIML class	E1
Permissible tolerance according to OIML (+/-)	0,008 mg

## Construction

Shape weight	Wire
Dimensions weight (W×D×H)	15×15×4 mm
Material	stainless steel polished

## Services (optional)

Article number for DAkkS calibration	962-259
--------------------------------------	---------

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	15×15×4 mm
Net weight	0,0005 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

## Optional cases

Case - Plastic article number	347-009-400
Case - Aluminium article number	317-009-600
Case - Wood article number	338-090-200

## Pictograms

### STANDARD

