

KERN 313-060-600



Aluminium protective box for weight sets (Cylindrical/compact, 1 mg - 1 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Aluminium protective box

Measuring System

Case suitable for OIML classes	E1 - M1
Suitable weight shape for case	Cylindrical/compact
Case suitable for nominal weight values	1 mg - 1 kg

Construction

Dimension (W×D×H)	230×140×110 mm
Material	aluminium
Case - Smallest possible weight	1 mg
Case - Highest possible weight	1 kg

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	230×140×110 mm
Shipping method	Parcel service
Net weight approx.	0,80 kg
Gross weight approx.	0,80 kg
Shipping weight	0,8 kg

Product information

GTIN/EAN number	4045761222793
REACH classification	no prohibited substances

Pictograms

STANDARD

