

KERN 313-07

KERN

Weight set, E2, 1 mg - 2 kg, Knob, stainless steel polished, in wooden case



Optional cases

Case - Wood article number 313-070-100

Product information

GTIN/EAN number 4045761003187
REACH classification no prohibited substances

Pictograms

STANDARD



OPTION



Category

| | |
|------------------|------------|
| Brand | KERN |
| Product category | Weight set |
| Product group | OIML E2 |
| Product family | 313 |

Measuring System

| | |
|-----------------------------|-------------|
| Weight OIML class | E2 |
| Test weight - nominal value | 1 mg - 2 kg |
| OIML-compliant | ✓ |

Construction

| | |
|--------------------|--------------------------|
| Dimension (W×D×H) | 305×190×140 mm |
| Shape weight | Knob |
| Material | stainless steel polished |
| Number of weights | 27 |
| Case of weight set | in wooden case |
| Adjustment method | Polishing |

Services (optional)

| | |
|--|---------|
| Article number for DAkkS calibration | 962-306 |
| Article number for guaranteed plus-deviation | 921-306 |

Packing & Shipping

| | |
|------------------------------|----------------|
| Delivery time | 1 d |
| Dimensions packaging (W×D×H) | 305×190×140 mm |
| Shipping method | Parcel service |
| Net weight approx. | 10 kg |
| Gross weight approx. | 10 kg |
| Shipping weight | 9,01 kg |