

# KERN 316-07

Individual weight, OIML Class E2, 100 g, Compact form, stainless steel polished



Article number for guaranteed plus-deviation 921-337

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	22×22×35 mm
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg

## Optional cases

Case - Plastic article number	317-070-400
Case - Aluminium article number	317-070-600
Case - Wood article number	317-070-100

## Product information

GTIN/EAN number	4045761094895
REACH classification	no prohibited substances

## Pictograms

### STANDARD



### OPTION



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	316

## Measuring System

Test weight - volume (cm <sup>3</sup> )	12,5
Conformity assessment possible	no
Test weight - Density uncertainty (kg/m <sup>3</sup> )	140
Test weight - volume uncertainty (cm <sup>3</sup> )	0,1563
Weight OIML class	E2
Test weight - nominal value	100 g
Permissible tolerance according to OIML (+/-)	0,16 mg
Test weight - Density (kg/m <sup>3</sup> )	7950
OIML-compliant	✓
Test weight - Density [Min] (OIML)	7810
Test weight - Density [Max] (OIML)	8210

## Construction

Shape weight	Compact form
Dimensions (Ø×H)	22×35 mm
Material	stainless steel polished
Adjustment method	Polishing

## Services (optional)

Article number for DAkkS calibration	962-337
--------------------------------------	---------