

KERN 317-05

Individual weight, OIML Class E2, 20 g, Knob, stainless steel polished



Article number for guaranteed plus-deviation	921-335
--	---------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	13×13×22 mm
Net weight	0,02 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	317-050-400
Case - Aluminium article number	317-050-600
Case - Wood article number	317-050-100

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	317

Measuring System

Conformity assessment possible	✓
Test weight - volume (cm ³)	2,5
Test weight - volume uncertainty (cm ³)	0,0313
Test weight - Density (kg/m ³)	7950
Test weight - Density uncertainty (kg/m ³)	140
Test weight - Density [Min] (OIML)	7500
Test weight - Density [Max] (OIML)	8570
Test weight - nominal value	20 g
OIML-compliant	✓
Weight OIML class	E2
Permissible tolerance according to OIML (+/-)	0,08 mg

Construction

Shape weight	Knob
Dimensions (Ø×H)	13×22 mm
Material	stainless steel polished

Services (optional)

Article number for conformity assessment	952-335
Article number for DAkkS calibration	962-335