

KERN 318-08

Milligram weights, E2, 200 mg, Platelet, stainless steel polished



Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	9×9×4 mm
Net weight	0,0002 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	347-009-400
Case - Aluminium article number	317-009-600
Case - Wood article number	338-090-200

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	318

Measuring System

Conformity assessment possible	✓
Test weight - volume (cm ³)	0,0233
Test weight - volume uncertainty (cm ³)	0,0005
Test weight - Density (kg/m ³)	7960
Test weight - Density uncertainty (kg/m ³)	170
Test weight - nominal value	200 mg
OIML-compliant	✓
Weight OIML class	E2
Permissible tolerance according to OIML (+/-)	0,02 mg

Construction

Shape weight	Platelet
Dimensions weight (W×D×H)	9×9×4 mm
Material	stainless steel polished

Services (optional)

Article number for conformity assessment	952-358
Article number for DAkkS calibration	962-358
Article number for guaranteed plus-deviation	921-358