

KERN 318-226

KERN

Set of milligram weights, E2, 1 mg - 500 mg, Platelet, stainless steel polished, in aluminium case



Gross weight approx.	0,40 kg
Shipping weight	0,39 kg

Optional cases

Case - Aluminium article number	313-010-600
---------------------------------	-------------

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML E2
Product family	318

Measuring System

Test weight - nominal value	1 mg - 500 mg
OIML-compliant	✓
Weight OIML class	E2

Construction

Dimension (W×D×H)	165×158×37 mm
Shape weight	Platelet
Material	stainless steel polished
Number of weights	13
Case of weight set	in aluminium case

Services (optional)

Article number for conformity assessment	952-350
Article number for DAkkS calibration	962-350
Article number for guaranteed plus-deviation	921-350

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	165×158×37 mm
Net weight	0,31 kg
Shipping method	Parcel service
Net weight approx.	0,35 kg