

# KERN 323-064

Weight set, F1, 1 mg - 1 kg, Knob, stainless steel polished, in plastic case



Gross weight approx.	3,8 kg
Shipping weight	4 kg

## Optional cases

Case - Plastic article number	313-082-400
-------------------------------	-------------

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	323

## Measuring System

Conformity assessment possible	✓
Test weight - nominal value	1 mg - 1 kg
OIML-compliant	✓
Weight OIML class	F1

## Construction

Dimension (W×D×H)	350×320×180 mm
Shape weight	Knob
Material	stainless steel polished
Number of weights	25
Case of weight set	in plastic case

## Services (optional)

Article number for conformity assessment	952-405
Article number for DAkkS calibration	962-405
Article number for guaranteed plus-deviation	921-405

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	350×320×180 mm
Shipping method	Parcel service
Net weight approx.	3,8 kg