

KERN 325-086

KERN

Weight set, F1, 1 mg - 5 kg, Eco Shape, stainless steel polished, in aluminium case



Optional cases

Case - Aluminium article number 313-080-600

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	325

Measuring System

Test weight - nominal value	1 mg - 5 kg
OIML-compliant	✓
Weight OIML class	F1

Construction

Dimension (W×D×H)	460×260×200 mm
Shape weight	Eco Shape
Material	stainless steel polished
Case of weight set	in aluminium case

Services (optional)

Article number for DAkkS calibration	962-407
Article number for guaranteed plus-deviation	921-407

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	460×260×200 mm
Net weight	17 kg
Shipping method	Parcel service
Net weight approx.	17 kg
Gross weight approx.	18 kg
Shipping weight	18 kg