

KERN 326-056

KERN

Weight set, F1, 1 g - 500 g, Eco Shape, stainless steel polished, in aluminium case



Optional cases

Case - Aluminium article number 314-050-600

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	326

Measuring System

Test weight - nominal value	1 g - 500 g
OIML-compliant	✓
Weight OIML class	F1

Construction

Dimension (W×D×H)	210×140×98 mm
Shape weight	Eco Shape
Material	stainless steel polished
Case of weight set	in aluminium case

Services (optional)

Article number for DAkkS calibration	962-418
Article number for guaranteed plus-deviation	921-418

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	210×140×98 mm
Net weight	0,71 kg
Shipping method	Parcel service
Net weight approx.	0,75 kg
Gross weight approx.	2,0 kg
Shipping weight	2 kg