

# KERN 326-36

# KERN

Individual weight, OIML Class F1, 5 kg, Block, stainless steel polished (OIML)



## Services (optional)

Article number for DAkkS calibration	962-443
--------------------------------------	---------

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	326

## Measuring System

Conformity assessment possible	✓
Test weight - nominal value	5 kg
Test weight - volume uncertainty (cm <sup>3</sup> )	7,8125
Permissible tolerance according to OIML (+/-)	25 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

## Construction

Shape weight	Block
Dimensions weight (W×D×H)	150×68×96 mm
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	150×68×96 mm
Net weight	5 kg
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	5,0 kg
Shipping weight	5 kg