

KERN 327-13

Individual weight, OIML Class F1, 5 kg, Knob, stainless steel polished



Article number for guaranteed plus-deviation	921-443
----------------------------------------------	---------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	80×80×147 mm
Net weight	5 kg
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	5,0 kg
Shipping weight	5 kg

Optional cases

Case - Plastic article number	347-130-400
Case - Aluminium article number	317-130-600
Case - Wood article number	317-130-100

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	327

Measuring System

Conformity assessment possible	✓
Test weight - volume (cm ³)	625
Test weight - volume uncertainty (cm ³)	7,8125
Test weight - Density (kg/m ³)	7950
Test weight - Density uncertainty (kg/m ³)	140
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Test weight - nominal value	5 kg
OIML-compliant	✓
Weight OIML class	F1
Permissible tolerance according to OIML (+/-)	25 mg

Construction

Shape weight	Knob
Dimensions (Ø×H)	80×147 mm
Material	stainless steel polished

Services (optional)

Article number for conformity assessment	952-443
Article number for DAkkS calibration	962-443