

KERN 327-16

KERN

Individual weight, OIML Class F1, 50 kg, Knob, stainless steel polished

Case - Wood article number 317-160-100

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	327

Measuring System

Conformity assessment possible	✓
Test weight - volume (cm ³)	6250
Test weight - volume uncertainty (cm ³)	78,125
Test weight - Density (kg/m ³)	7950
Test weight - Density uncertainty (kg/m ³)	140
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Test weight - nominal value	50 kg
OIML-compliant	✓
Weight OIML class	F1
Permissible tolerance according to OIML (+/-)	250 mg

Construction

Shape weight	Knob
Dimensions (Ø×H)	178×298 mm
Material	stainless steel polished

Services (optional)

Article number for conformity assessment	952-446
Article number for DAkkS calibration	962-446
Article number for guaranteed plus-deviation	921-446

Packing & Shipping

Delivery time	2 d
Dimensions packaging (W×D×H)	178×178×298 mm
Net weight	50 kg
Shipping method	Freight
Net weight approx.	50 kg
Gross weight approx.	50 kg
Shipping weight	50 kg

Optional cases

Case - Aluminium article number	317-160-600
---------------------------------	-------------