


KERN 328-01

KERN

Milligram weights, F1, 1 mg, Platelet, stainless steel polished



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	328

Measuring System

Conformity assessment possible	✓
Test weight - volume (cm ³)	0,0004
Test weight - volume uncertainty (cm ³)	0,0001
Test weight - Density (kg/m ³)	7960
Test weight - Density uncertainty (kg/m ³)	170
Test weight - nominal value	1 mg
OIML-compliant	✓
Weight OIML class	F1
Permissible tolerance according to OIML (+/-)	0,02 mg

Construction

Shape weight	Platelet
Dimensions weight (W×D×H)	3,5×3,5×3 mm
Material	stainless steel polished

Services (optional)

Article number for conformity assessment	952-451
Article number for DAkkS calibration	962-451
Article number for guaranteed plus-deviation	921-451

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	3,5×3,5×3 mm
Net weight	0,01 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	347-009-400
Case - Aluminium article number	317-009-600
Case - Wood article number	338-090-200

Pictograms

STANDARD

