

KERN 328-38

KERN

Milligram weights, F1, 200 mg, Wire, stainless steel polished



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	328

Measuring System

Conformity assessment possible	✓
Test weight - volume (cm ³)	0,025
Test weight - volume uncertainty (cm ³)	0,0005
Test weight - Density (kg/m ³)	8000
Test weight - Density uncertainty (kg/m ³)	160
Test weight - nominal value	200 mg
OIML-compliant	✓
Weight OIML class	F1
Permissible tolerance according to OIML (+/-)	0,06 mg

Construction

Shape weight	Wire
Dimensions weight (W×D×H)	10×10×4 mm
Material	stainless steel polished

Services (optional)

Article number for conformity assessment	952-458
Article number for DAkkS calibration	962-458
Article number for guaranteed plus-deviation	921-458

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	10×10×4 mm
Net weight	0,0002 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	347-009-400
Case - Aluminium article number	317-009-600
Case - Wood article number	338-090-200

Pictograms

STANDARD

