

# KERN 334-026

# KERN

Weight set, F2, 1 g - 50 g, Knob, stainless steelfine turned, in aluminium case



## Optional cases

Case - Aluminium article number 314-020-600

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	334

## Measuring System

Test weight - nominal value	1 g - 50 g
OIML-compliant	✓
Weight OIML class	F2

## Construction

Dimension (W×D×H)	165×120×75 mm
Shape weight	Knob
Material	stainless steelfine turned
Case of weight set	in aluminium case

## Services (optional)

Article number for DAkkS calibration	962-415
Article number for guaranteed plus-deviation	921-415

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	165×120×75 mm
Net weight	0,8 kg
Shipping method	Parcel service
Net weight approx.	0,80 kg
Gross weight approx.	0,80 kg
Shipping weight	0,8 kg