

KERN 334-036

KERN

Weight set, F2, 1 g - 100 g, Knob, stainless steelfine turned, in aluminium case



Optional cases

Case - Aluminium article number 314-030-600

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	334

Measuring System

Test weight - nominal value	1 g - 100 g
OIML-compliant	✓
Weight OIML class	F2

Construction

Dimension (W×D×H)	210×140×98 mm
Shape weight	Knob
Material	stainless steelfine turned
Case of weight set	in aluminium case

Services (optional)

Article number for DAkkS calibration	962-416
Article number for guaranteed plus-deviation	921-416

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	210×140×98 mm
Net weight	0,9 kg
Shipping method	Parcel service
Net weight approx.	0,90 kg
Gross weight approx.	0,90 kg
Shipping weight	0,9 kg