

KERN 337-02

KERN

Individual weight, OIML Class F2, 2 g, Knob, stainless steelfine turned



Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	6×6×11 mm
Net weight	0,002 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	347-030-400
Case - Aluminium article number	317-020-600
Case - Wood article number	337-020-200

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F2
Product family	337

Measuring System

Conformity assessment possible	no
Test weight - volume (cm ³)	0,25
Test weight - volume uncertainty (cm ³)	0,0031
Test weight - Density (kg/m ³)	7950
Test weight - Density uncertainty (kg/m ³)	140
Test weight - Density [Min] (OIML)	2000
Test weight - nominal value	2 g
OIML-compliant	✓
Weight OIML class	F2
Permissible tolerance according to OIML (+/-)	0,4 mg

Construction

Shape weight	Knob
Dimensions (Ø×H)	6×11 mm
Material	stainless steelfine turned

Services (optional)

Article number for DAkkS calibration	962-432
Article number for guaranteed plus-deviation	921-432