

KERN 337-11



Individual weight, OIML Class F2, 1 kg, Knob, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration	962-441
Article number for guaranteed plus-deviation	921-441

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F2
Product family	337

Measuring System

Conformity assessment possible	no
Test weight - nominal value	1 kg
Test weight - volume (cm ³)	125
Test weight - volume uncertainty (cm ³)	1,5625
Permissible tolerance according to OIML (+/-)	16 mg
Test weight - Density [Min] (OIML)	6400
Test weight - Density [Max] (OIML)	10700
Weight OIML class	F2

Construction

Shape weight	Knob
Dimensions (Ø×H)	48×80 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	48×48×80 mm
Net weight	1 kg
Shipping method	Parcel service
Net weight approx.	1,00 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg