

KERN 337-16

Individual weight, OIML Class F2, 50 kg, Knob, stainless steel fine turned



Article number for guaranteed plus-deviation 921-446

Packing & Shipping

Delivery time	17 d
Dimensions packaging (W×D×H)	178×178×298 mm
Shipping method	Freight
Net weight approx.	50 kg
Gross weight approx.	50 kg
Shipping weight	50 kg

Optional cases

Case - Aluminium article number	317-160-600
Case - Wood article number	337-160-200

Product information

GTIN/EAN number	404576 1006638
REACH classification	no prohibited substances

Pictograms

STANDARD



OPTION



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F2
Product family	337

Measuring System

Test weight - volume (cm ³)	6250
Conformity assessment possible	no
Test weight - Density uncertainty (kg/m ³)	140
Test weight - volume uncertainty (cm ³)	78,125
Weight OIML class	F2
Test weight - nominal value	50 kg
Permissible tolerance according to OIML (+/-)	800 mg
Test weight - Density (kg/m ³)	7950
OIML-compliant	✓
Test weight - Density [Min] (OIML)	6400
Test weight - Density [Max] (OIML)	10700

Construction

Shape weight	Knob
Dimensions (Ø×H)	178×298 mm
Material	stainless steel fine turned
Adjustment method	Adjusting cavity - Seal

Services (optional)

Article number for DAkkS calibration	962-446
--------------------------------------	---------