

# KERN 343-02

Weight set, M1, 1 mg - 50 g, Knob, stainless steel fine turned, in wooden case



## Optional cases

Case - Wood article number 353-420-200

## Product information

GTIN/EAN number 4045761007468  
REACH classification no prohibited substances

## Pictograms

### STANDARD



### OPTION



## Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	343

## Measuring System

Weight OIML class	M1
Test weight - nominal value	1 mg - 50 g
OIML-compliant	✓

## Construction

Dimension (W×D×H)	129×86×40 mm
Shape weight	Knob
Material	stainless steel fine turned
Number of weights	18
Case of weight set	in wooden case

## Services (optional)

Article number for DAkkS calibration	962-601
--------------------------------------	---------

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	129×86×40 mm
Shipping method	Parcel service
Net weight approx.	0,40 kg
Gross weight approx.	0,40 kg
Shipping weight	0,38 kg