

KERN 343-066

Weight set, M1, 1 mg - 1 kg, Knob, stainless steelfine turned, in aluminium case



Optional cases

Case - Aluminium article number 313-060-600

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	343

Measuring System

Test weight - nominal value	1 mg - 1 kg
OIML-compliant	✓
Weight OIML class	M1

Construction

Dimension (W×D×H)	230×140×110 mm
Shape weight	Knob
Material	stainless steelfine turned
Case of weight set	in aluminium case

Services (optional)

Article number for conformity assessment	952-605
Article number for DAkkS calibration	962-605

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	230×140×110 mm
Net weight	2,71 kg
Shipping method	Parcel service
Net weight approx.	2,8 kg
Gross weight approx.	3,6 kg
Shipping weight	3,5 kg