

KERN 343-06

KERN

Weight set, M1, 1 mg - 1 kg, Knob, stainless steelfine turned, in wooden case



Optional cases

Case - Wood article number 353-460-200

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	343

Measuring System

Test weight - nominal value	1 mg - 1 kg
OIML-compliant	✓
Weight OIML class	M1

Construction

Dimension (W×D×H)	230×180×112 mm
Shape weight	Knob
Material	stainless steelfine turned
Number of weights	25
Case of weight set	in wooden case

Services (optional)

Article number for conformity assessment	952-605
Article number for DAkkS calibration	962-605

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	230×180×112 mm
Net weight	4,21 kg
Shipping method	Parcel service
Net weight approx.	4,4 kg
Gross weight approx.	4,4 kg
Shipping weight	4,21 kg