

KERN 347-015

KERN

Individual weight, OIML Class M1, 1 g, Slot weight (Type A), stainless steelfine turned



Net weight	0,001 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	347-030-400
Case - Aluminium article number	nicht verfügbar

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System

Test weight - volume (cm ³)	0,119
Test weight - volume uncertainty (cm ³)	0,0024
Test weight - Density (kg/m ³)	8000
Test weight - Density uncertainty (kg/m ³)	100
Test weight - Density [Min] (OIML)	2600
Test weight - nominal value	1 g
OIML-compliant	✓
Weight OIML class	M1
Permissible tolerance according to OIML (+/-)	1 mg

Construction

Shape weight	Slot weight (Type A)
Dimensions (Ø×H)	13×2 mm
Material	stainless steelfine turned

Services (optional)

Article number for DAkkS calibration	962-631
--------------------------------------	---------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	13×13×2 mm