

KERN 347-03



Individual weight, OIML Class M1, 5 g, Knob, stainless steelfine turned (OIML)



Services (optional)

Article number for conformity assessment	952-633
Article number for DAkkS calibration	962-633

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System

Conformity assessment possible	✓
Test weight - nominal value	5 g
Test weight - volume (cm ³)	0,625
Test weight - volume uncertainty (cm ³)	0,0078
Permissible tolerance according to OIML (+/-)	1,6 mg
Weight OIML class	M1

Construction

Shape weight	Knob
Dimensions (Ø×H)	8×16,5 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	8×8×16 mm
Net weight	0,01 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg