

KERN 347-05



Individual weight, OIML Class M1, 20 g, Knob, stainless steelfine turned (OIML)



Services (optional)

Article number for conformity assessment	952-635
Article number for DAkkS calibration	962-635

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System

Conformity assessment possible	✓
Test weight - nominal value	20 g
Test weight - volume (cm ³)	2,5
Test weight - volume uncertainty (cm ³)	0,0313
Permissible tolerance according to OIML (+/-)	2,5 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

Construction

Shape weight	Knob
Dimensions (Ø×H)	13×23,5 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	13×13×23 mm
Net weight	0,02 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg