

KERN 347-065



Individual weight, OIML Class M1, 50 g, Slot weight, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration 962-636

Pictograms

Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System

Test weight - nominal value	50 g
Test weight - volume (cm ³)	5,9524
Test weight - volume uncertainty (cm ³)	0,1205
Permissible tolerance according to OIML (+/-)	3 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

Construction

Shape weight	Slot weight
Dimensions (Ø×H)	28×12 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	28×28×12 mm
Net weight	0,05 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg