

KERN 347-125

KERN

Individual weight, OIML Class M1, 2 kg, Slot weight (Type A), stainless steelfine turned



Net weight	2 kg
Shipping method	Parcel service
Net weight approx.	2,0 kg
Gross weight approx.	2,0 kg
Shipping weight	2 kg

Optional cases

Case - Plastic article number	347-130-400
-------------------------------	-------------

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

Measuring System

Test weight - volume (cm ³)	238,0952
Test weight - volume uncertainty (cm ³)	4,8186
Test weight - Density (kg/m ³)	8000
Test weight - Density uncertainty (kg/m ³)	100
Test weight - Density [Min] (OIML)	2600
Test weight - nominal value	2 kg
OIML-compliant	✓
Weight OIML class	M1
Permissible tolerance according to OIML (+/-)	100 mg

Construction

Shape weight	Slot weight (Type A)
Dimensions (Ø×H)	80×57 mm
Material	stainless steelfine turned

Services (optional)

Article number for DAkkS calibration	962-642
--------------------------------------	---------

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	80×80×57 mm