

# KERN 347-146



Individual weight, OIML Class M1, 10 kg, Hook weight, stainless steelfine turned (OIML)



## Services (optional)

Article number for DAkkS calibration	962-644
--------------------------------------	---------

## Pictograms

### Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	347

### Measuring System

Test weight - nominal value	10 kg
Test weight - volume (cm <sup>3</sup> )	1190,4762
Test weight - volume uncertainty (cm <sup>3</sup> )	24,093
Permissible tolerance according to OIML (+/-)	500 mg
Test weight - Density [Min] (OIML)	2600
Weight OIML class	M1

### Construction

Shape weight	Hook weight
Dimensions (Ø×H)	220×160 mm
Material	stainless steelfine turned (OIML)

### Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	220×220×160 mm
Net weight	10 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	10 kg