

KERN 348-08

KERN

Milligram weights, M1, 200 mg, Platelet, stainless steel polished



Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg

Optional cases

Case - Plastic article number	347-009-400
Case - Aluminium article number	317-009-600
Case - Wood article number	338-090-200

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	348

Measuring System

Conformity assessment possible	✓
Test weight - Density (kg/m ³)	7960
Test weight - Density uncertainty (kg/m ³)	170
Test weight - nominal value	200 mg
OIML-compliant	✓
Weight OIML class	M1
Permissible tolerance according to OIML (+/-)	0,6 mg

Construction

Shape weight	Platelet
Dimensions weight (W×D×H)	9×9×4 mm
Material	stainless steel polished

Services (optional)

Article number for conformity assessment	952-658
Article number for DAkkS calibration	962-658

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	9×9×4 mm
Net weight	0,0002 kg