

# KERN 353-490-200



F2 + M1 + M2 Wooden box for individual weight sets 1 mg - 10 kg



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

## Measuring System

Case suitable for OIML classes	F2 + M1 + M2
Case suitable for nominal weight values	1 mg - 10 kg
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (W×D×H)	540×190×245 mm
Material	beech

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	540×190×245 mm
Net weight	8,8 kg
Shipping method	Parcel service
Net weight approx.	9 kg
Gross weight approx.	11 kg
Shipping weight	10,5 kg

## Pictograms

STANDARD

