

# KERN 366-94

**KERN**

Individual weight, OIML Class M3, 1 kg, Knob, cast iron lacquered



Dimensions packaging (W×D×H)	55×55×87 mm
Net weight	1 kg
Shipping method	Parcel service
Net weight approx.	1,00 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	366

## Measuring System

Conformity assessment possible	no
Test weight - volume (cm <sup>3</sup> )	140,8451
Test weight - volume uncertainty (cm <sup>3</sup> )	11,9024
Test weight - Density (kg/m <sup>3</sup> )	7100
Test weight - Density uncertainty (kg/m <sup>3</sup> )	600
Test weight - Density [Min] (OIML)	1500
Test weight - nominal value	1 kg
OIML-compliant	✓
Weight OIML class	M3
Permissible tolerance according to OIML (+/-)	500 mg

## Construction

Shape weight	Knob
Dimensions (Ø×H)	55×87 mm
Material	cast iron lacquered

## Services (optional)

Article number for DAkkS calibration	962-641
--------------------------------------	---------

## Packing & Shipping

Delivery time	1 d
---------------	-----