

# KERN 367-02



Individual weight, OIML Class M3, 2 g, Knob, stainless steelfine turned (OIML)



## Services (optional)

Article number for DAkkS calibration	962-632
--------------------------------------	---------

## Pictograms

## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	367

## Measuring System

Conformity assessment possible	no
Test weight - nominal value	2 g
Test weight - volume (cm <sup>3</sup> )	0,2381
Test weight - volume uncertainty (cm <sup>3</sup> )	0,0048
Permissible tolerance according to OIML (+/-)	12 mg
Weight OIML class	M3

## Construction

Shape weight	Knob
Dimensions (Ø×H)	6×11 mm
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	6×6×11 mm
Net weight	0,002 kg
Shipping method	Parcel service
Net weight approx.	0,05 kg
Gross weight approx.	0,05 kg
Shipping weight	0,1 kg