

# KERN 367-11



Individual weight, OIML Class M3, 1 kg, Knob, stainless steelfine turned (OIML)



## Services (optional)

Article number for DAkkS calibration	962-641
--------------------------------------	---------

## Pictograms

## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	367

## Measuring System

Conformity assessment possible	no
Test weight - nominal value	1 kg
Test weight - volume (cm <sup>3</sup> )	119,0476
Test weight - volume uncertainty (cm <sup>3</sup> )	2,4093
Permissible tolerance according to OIML (+/-)	500 mg
Test weight - Density [Min] (OIML)	1500
Weight OIML class	M3

## Construction

Shape weight	Knob
Dimensions (Ø×H)	48×80 mm
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	48×48×80 mm
Net weight	1 kg
Shipping method	Parcel service
Net weight approx.	1,00 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg