

# KERN 367-12



Individual weight, OIML Class M3, 2 kg, Knob, stainless steelfine turned (OIML)



## Services (optional)

Article number for DAkkS calibration	962-642
--------------------------------------	---------

## Pictograms

### Category

Brand	KERN
Product category	Individual weight
Product group	OIML M3
Product family	367

### Measuring System

Conformity assessment possible	no
Test weight - nominal value	2 kg
Test weight - volume (cm <sup>3</sup> )	238,0952
Test weight - volume uncertainty (cm <sup>3</sup> )	4,8186
Permissible tolerance according to OIML (+/-)	1 g
Test weight - Density [Min] (OIML)	1500
Weight OIML class	M3

### Construction

Shape weight	Knob
Dimensions (Ø×H)	60×100 mm
Material	stainless steelfine turned (OIML)

### Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	60×60×100 mm
Net weight	2 kg
Shipping method	Parcel service
Net weight approx.	2,0 kg
Gross weight approx.	2,0 kg
Shipping weight	2 kg