

# KERN MAP 80K1S

# KERN

Hand grip dynamometer, e.g. for rehabilitation treatment after accidents



## Category

Brand	KERN
Product category	Medical scale
Product group	Hand grip dynamometer
Product family	MAP

## Measuring System

Weighing capacity [Max]	80 kg
Readability [d]	100 g
Recommended adjusting weight	40 kg (M1)
Possible calibration points	40 kg
Units	kg lb

## Display

Display type	LCD graphic
Display digit height	12 mm

## Construction

Dimensions device (W×D×H)	55×88×212 mm
---------------------------	--------------

## Functions

Auto-off interval(s) in battery mode/rechargeable battery mode	1 min
--	-------

## Power Supply

Supplied power supply	Battery
Battery / accumulator type	Li-Metal
Battery	CR2450
Battery operating time	53 h
Battery capacity	630 mAh
Battery voltage	3 V

## Environmental conditions

Humidity of environment [Max]	80 %
Use temperature [Min]	5 °C
Use temperature [Max]	35 °C
Storage temperature [Min]	-20 °C
Storage temperature [Max]	60 °C

## Services (optional)

Article number for factory calibration	961-102K
--	----------

## Packing & Shipping

Dimensions transport container (W×D×H)	350×265×85 mm
Delivery time	1 d
Dimensions packaging (W×D×H)	367×285×100 mm
Net weight	0,31 kg
Shipping method	Parcel service
Net weight approx.	0,35 kg
Gross weight approx.	1,8 kg
Shipping weight	2,1 kg

## Pictograms

### STANDARD



### OPTION

