

KERN MFB 150K100S05

KERN

Design body analysis scale to determine body weight, body fat, body fluid, muscle mass, proportion of bone, BMI and metabolic rate/metabolic rate when exercising



Dimensions weighing surface (W×D)	300×300 mm
Dimensions weighing surface	300×300 mm
Dimensions weighing platform (W×D×H)	300×300×25 mm
Material weighing plate	glass

Functions

Auto-off interval(s) in battery mode/rechargeable battery mode	1 min
--	-------

Power Supply

Supplied power supply	Battery
Power supply cable removable	removable
Battery / accumulator type	Li-Metal
Battery	2×CR2032
Battery operating time	100 h
Battery capacity	210 mAh
Battery – Number of cells	10
Battery – Individual weight	30 g
Battery – Total weight	300 g
Battery – Construction	Individual cells
Battery capacity (Wh)	0,63 Wh
Battery – Lithium content	550 mg
Battery – Tests passed UN 38.3	✓
Battery/Accu connected	not connected, isolated in the device

Environmental conditions

Ambient temperature [Min]	-10 °C
Ambient temperature [Max]	40 °C
Humidity of environment [Min]	5 %
Humidity of environment [Max]	80 %
Storage temperature [Min]	-20 °C
Storage temperature [Max]	60 °C

Approval

CE mark	✓
---------	---

Services (optional)

Article number for DAkkS calibration	963-129S05
Article number for certificate of conformity	969-517

Packing & Shipping

Delivery time	1 d
---------------	-----

Category

Brand	KERN
Product category	Medical scale
Product group	Body analysis scale
Product family	MFB

Measuring System

Weighing system	Strain gauge
Weighing capacity [Max]	182 kg
Readability [d]	0,1 kg
Reproducibility	0,1 kg
Linearity	± 0,5 kg
Resolution	1.820
Adjustment options	Calibration not possible
Stabilization time	2 s
Warm-up time	10 min
Eccentric loading at 1/3 [Max]	0,5 kg
Construction type of the scale	Single-range balance
Units	kg lb
Maximum creep (15 minutes)	200 g
Maximum creep (30 minutes)	400 g
Default unit	kg

Display

Display type	LCD
Display backlight	✓
Display digit height	35 mm

Construction

Dimensions completely mounted (W×D×H)	300×300×25 mm
---------------------------------------	---------------

KERN MFB 150K100S05



Design body analysis scale to determine body weight, body fat, body fluid, muscle mass, proportion of bone, BMI and metabolic rate/metabolic rate when exercising

Dimensions packaging (W×D×H) 330×330×300 mm

Shipping method Parcel service

Net weight approx. 9 kg

Gross weight approx. 10 kg

Shipping weight 10 kg

Product information

GTIN/EAN number 404576 1208360

REACH classification no prohibited substances

Pictograms

STANDARD

