


KERN MPK 200K-1P

KERN

Personal scale with column for convenient, ergonomically-optimised operation



Category

Brand	KERN
Product category	Medical scale
Product group	Personal floor scale
Product family	MPK

Measuring System

Construction type of the scale	Single-range balance
Weighing system	Strain gauge
Weighing capacity [Max]	250 kg
Readability [d]	0,1 kg
Reproducibility	0,1 kg
Linearity	± 1 kg
Resolution	2.500
Adjustment options	Adjusting with external weight
Recommended adjusting weight	250 kg (M1)
Stabilization time	2 s
Warm-up time	10 min
Eccentric loading at 1/3 [Max]	0,5 kg
Default unit	kg
Units	kg lb

Display

Display type	LCD
Display digit height	28 mm
Languages of the user interface	Symbol language

Construction

Dimension housing (W×D×H)	350×370×960 mm
---------------------------	----------------

Dimensions completely mounted (W×D×H)	350×420×955 mm
Dimension display device (W×D×H)	190×73×27,5 mm
Dimensions weighing surface (W×D)	350×350 mm
Dimensions weighing surface	350×350 mm
Dimensions weighing platform (W×D×H)	350×350×45 mm
Material housing	ABS plastic
Material weighing plate	plastic
Material display housing	Plastics
Stand height	900 mm

Functions

Hold function	✓
Auto-off interval(s) in battery mode/rechargeable battery mode	10 sec

Power Supply

Supplied power supply	Battery
Battery / accumulator type	Alkali(-Manganese)
Battery	4×1.5 V AA
Battery connection	Battery insert
Battery operating time	100 h
Battery voltage	6 V

Environmental conditions

Ambient temperature [Min]	10 °C
Ambient temperature [Max]	40 °C
Humidity of environment [Min]	5 %
Humidity of environment [Max]	80 %

Approval

CE mark	✓
---------	---

Services (optional)

Article number for adjustment at the location of installation	961-249
Article number for DAkkS calibration	963-129
Article number for certificate of conformity (verification)	969-517

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	445×1003×145 mm
Net weight	3,8 kg
Shipping method	Parcel service

KERN MPK 200K-1P



Personal scale with column for convenient, ergonomically-optimised operation

Net weight approx.	3,8 kg
Gross weight approx.	7 kg
Shipping weight	12,9 kg

Pictograms

STANDARD



OPTION

