

# KERN OBB-A 1441

**KERN**

Achromatic objective lens, 100 x, Plan / 1 (water) (sprung) W.D. (0,18 mm)



## Pictograms

### STANDARD



## Category

|                  |                                 |
|------------------|---------------------------------|
| Brand            | Optics                          |
| Product category | Optical Instrument              |
| Product group    | Microscope objective achromatic |

## Construction

|             |   |
|-------------|---|
| Anti fungus | ✓ |
|-------------|---|

## Objective

|                            |                 |
|----------------------------|-----------------|
| Objective magnification    | 100 x           |
| Objective quality          | Plan achromatic |
| Numerical Aperture         | 1               |
| Objective working distance | 0,18 mm         |

## Environmental conditions

|                           |       |
|---------------------------|-------|
| Storage temperature [Min] | -5 °C |
| Storage temperature [Max] | 40 °C |

## Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Delivery time                | 1 d            |
| Dimensions packaging (W×D×H) | 65×40×40 mm    |
| Net weight                   | 0,3 kg         |
| Shipping method              | Parcel service |
| Net weight approx.           | 0,30 kg        |
| Gross weight approx.         | 0,60 kg        |
| Shipping weight              | 0,6 kg         |