

KERN OBB-A 1632



3W LED Fluorescence unit (B / G / V / UV) for inverted microscope



Category

Brand	Optics
Product category	Optical Instrument
Product group	Fluorescence unit

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	410×390×310 mm
Net weight	8 kg
Shipping method	Parcel service
Net weight approx.	8 kg
Gross weight approx.	9 kg
Shipping weight	9,9 kg

Pictograms