

KERN OBS 114

KERN

The school microscope – for the first steps in microscopy and for use in biology lessons



Category

Brand	Optics
Product category	Microscope
Product group	Compound microscope
Product family	OBS-1

Construction

Dimension (W×D×H)	130×300×310 mm
Optical system	Finite
Tube type	Binocular
Tube angle of inclination	30°
Tube 360° rotation	✓
Tube type	Siedentopf
Contrasting methods	Bright field
	4×
Standard objectives	10×
	40×
Lens quality	achromatic
Nosepiece screw-in locations	4
Eyepieces fixed	✓
Interpupillary distance [Min]	55 mm
Interpupillary distance [Max]	75 mm
Diopter adjustment	both-sided
Diopter adjustment [Min]	-5
Diopter adjustment [Max]	5

Ocular

Ocular type	Eyepiece WF 10 x / Ø 18mm with anti-fungus
Ocular field width	WF
Eye point	Standard
Ocular magnifications	10 x

Ocular visual field	18 mm
Ocular diameter	23,2 mm

Objective

Objectives - Details	Objective Achromatic 10 x / 0,25 anti-fungus Objective Achromatic 4 x / 0,1 anti-fungus Objective Achromatic 40 x / 0,65 spring, anti-fungus
Objectives magnifications	4 x / 10 x / 40 x

Illumination

Illuminance	Transmitted light
Illumination type transmitted light	LED
Illumination intensity transmitted light	0,5 W
Filter possible	✓
Aperture diaphragm	✓

Focussing

Field of view [Min]	0,45 mm
Field of view [Max]	4,5 mm
Fine drive minimum	0,00025 mm
Focusing mechanism	coaxial coarse and fine drive

Power Supply

Supplied power supply	Power supply unit
Plug-in power supply type	Power adapter
Plug-in power supply / adapter for countries - included with the delivery	EURO
Plug-in power supply / adapter for countries - optional	EURO UK
Input voltage power supply / power [Max]	100 - 240 V
Input voltage device / power [Max]	3V, 500mA

Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

Approval

CE mark	✓
---------	---

KERN OBS 114

The school microscope – for the first steps in microscopy and for use in biology lessons

Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	410×280×210 mm
Net weight	3,138 kg
Shipping method	Parcel service
Net weight approx.	3,2 kg
Gross weight approx.	4,0 kg
Shipping weight	4,8 kg

Pictograms

STANDARD

