

# KERN OBT 105

# KERN

The modern compound microscope for teaching in your class room



Category	
Brand	Optics
Product category	Microscope
Product group	Compound microscope
Product family	OBT-1

Construction	
Dimension (WxDxH)	195x147x325 mm
Optical system	Finite
Tube type	Monocular
Tube angle of inclination	45°
Tube 360° rotation	✓
Contrasting methods	Bright field
Standard objectives	4x 10x 40x 100x
Lens quality	achromatic
Nosepiece screw-in locations	4
Eyepieces fixed	✓

Ocular	
Ocular type	Eyepiece WF 10 x / Ø 18mm
Ocular field width	WF
Eye point	Standard
Ocular magnifications	10 x
Ocular visual field	18 mm
Ocular diameter	23,2 mm

Objective	
Objectives - Details	Objective Achromatic 10 x / 0,25 Objective Achromatic 100 x / 1,25 oil, spring Objective Achromatic 4 x / 0,1 Objective Achromatic 40 x / 0,65 spring
Objectives magnifications	4 x / 10 x / 40 x / 100 x

Illumination	
Illuminance	Transmitted light
Illumination type transmitted light	LED
Illumination intensity transmitted light	1 W
Illumination dimmable	Transmitted light
Filter possible	✓
Aperture diaphragm	✓

Focussing	
Field of view [Min]	0,45 mm
Field of view [Max]	4,5 mm
Fine drive minimum	0,0025 mm
Focusing mechanism	coaxial coarse and fine drive

Power Supply	
Supplied power supply	Power supply & battery function (battery not included)
Plug-in power supply type	Power adapter
Plug-in power supply / adapter for countries - included with the delivery	EURO
Plug-in power supply / adapter for countries - optional	AUS UK US CH
Input voltage power supply / power [Max]	100 - 240 V
Input voltage device / power [Max]	5 V 1000 mA
Battery connection	Pad + Spring

Environmental conditions	
Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

# KERN OBT 105



The modern compound microscope for teaching in your class room

## Approval

CE mark ✓

## Packing & Shipping

Delivery time	1 d
Dimensions packaging (W×D×H)	400×280×195 mm
Net weight	2,755 kg
Shipping method	Parcel service
Net weight approx.	2,8 kg
Gross weight approx.	3,4 kg
Shipping weight	4,4 kg

## Pictograms

### STANDARD



### OPTION

