

# KERN OZB-A4604



Object clamp for stereo microscope stands



## Category

Brand	Optics
Product category	Optical Instrument
Product group	Object clamp

## Construction

Material	metal
----------	-------

## Environmental conditions

Storage temperature [Min]	-5 °C
Storage temperature [Max]	40 °C

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	130×45×15 mm
Net weight	0,1 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,30 kg
Shipping weight	0,3 kg

## Pictograms

STANDARD

